

## Free sciNote User Roles and Permissions

ACTIONS	TEAM			PROJECT				Comment
	Team admin	Normal user	Guest	Owner	User	Technician	Viewer	
<b>TEAM PERMISSIONS</b>								
Create projects	x	x						
Create protocols in repository	x	x						
Create repositories	x							
Create repository columns	x	x						
Create repository rows (items)	x	x						
Create sample columns	x	x						
Create samples	x	x						
Manage sample types and groups	x	x						Manage includes: create, edit, delete.
Manage team users	x							Manage includes: invite, edit role, remove user from the team.
Read team	x	x	x					
Update team	x							
<b>PROJECT PERMISSIONS</b>								
Archive project				x				
Create comments in project				x	x	x		
Create experiments				x	x			
Manage project				x				Manage includes: create, edit, archive.
Manage reports				x	x	x		Manage includes: create, edit, delete.
Manage tags				x	x			Manage includes: create, edit, (un)assign, delete.
Read project	x	x	x					You can see projects if you are a member of project. Project is public, or you are a team Administrator.
Restore project				x				
Manage comment in project				x				Manage includes: edit and delete. You can manage only your on comments, only Owner can manage all comments.

### EXPERIMENT PERMISSIONS

Archive experiment		x	x		
Clone experiment	x	x	x	x	
Manage experiment			x	x	Manage includes: create, edit.
Move experiment	x	x	x	x	
Read experiment					You can see experiments if you are a member of project. Project is public, or you are a team Administrator.
Restore experiment			x	x	

### TASK PERMISSIONS

Assign repository items to task			x	x	x	
Assign sample to task			x	x	x	
Complete task			x	x	x	
Create comments in task			x	x	x	
Manage task			x	x		Manage includes: create, edit, move, archive.
Manage users in task			x			Manage includes: (un)assign, edit role.
Restore task			x	x		
Manage comment in task			x			Manage includes: create, edit. You can manage only your on comments, only Owner can manage all comments.

### RESULT PERMISSIONS

Manage result			x			Manage includes: create, edit, archive, delete.
---------------	--	--	---	--	--	-------------------------------------------------

### PROTOCOL PERMISSIONS

Clone protocol in repository	x	x				
Complete step or check-box in a step			x	x	x	
Manage protocol in task			x	x		Manage includes: create, edit, delete, unlink, revert, update from repository, import.
Manage protocol in repository	x	x				Manage includes: create, edit, delete, reorder step, toggle visibility, and archive. You can only manage your own protocols.
Read protocol in task	x	x	x			
Read protocol in repository	x	x	x			
Restore protocol in repository	x	x				You can only manage your own protocols.

### SAMPLE PERMISSIONS

Manage sample	x	x				Manage includes: edit, delete.
Manage sample column	x	x				Manage includes: create, edit, delete.

### CUSTOM REPOSITORY PERMISSIONS

Manage repository	x		Manage includes: create, edit, delete.
Manage repository item	x	x	Manage includes: edit, delete.
Manage repository column	x	x	Manage includes: create, edit, delete.